

SWISSFLOAT

FOR ROUTINE EQUINE DENTISTRY

The #1 Motorized Float in Europe for over 20 Years

16" or 21" Slimline Model and 7" Incisor Model w/o Disk Guards

Ideal to reduce Hooks, Points and Incisors or Canines


Variable Speed 120V Corded or 12V/18V Cordless Motor

Fine Carbide and Durable Diamond Chip Grinding Disk

Gentle, Safe and Fast with Protecting Disk Guards


"The Work Horse in My Dentistry Practice"

"The Best Money Maker I've Ever Had"



SWISS QUALITY

SAFE GENTLE FAST



WWW.SWISSVET.COM

Distributed and Serviced in the US, Canada and Mexico

Swissfloat is designed for equine use only.

Although easy to handle Swissfloat must be used by qualified and horse experienced persons only. Always inform your client of possible adverse effects of dental procedures.

General precautions:

- If you have long hair, ensure that it is properly tied back.
- Wear protective gloves
- Ensure suitable eye protection.
- Consider ear protection as excessive noise may damage hearing
- Consider mask to protect you from potential airborne particles.
- Make sure that you are familiar with dentistry techniques
- Be always aware of the potential hazards of treatment procedures.
- 120V/220V apparatus should be used with GFCI protection only
- Keep materials, tools and hands dry while using electrical apparatus.
- Bulky clothing, rings, watches and dangling jewelry are dangerous when performing equine dentistry

Keep Swissfloat clean and dry after use. Lubrication is not required. Use water for cleaning restrictively.

No daily maintenance - but service is recommended once a year



Tel (877) 794 7735 - (877) swissfloat

FAX (865) 540 8850

mail@swissvet.com